

VT75 BIOMEDICAL CHEST FREEZER

The VT low temperature freezers creates the possibility to maintain temperatures as low as -60°C. Supreme stability, reliability, user-friendliness and ease of cleaning make these freezers an ideal solution for laboratories and hospitals.





| DIMENSIONS | |
|------------------------------------|----------------------------------|
| Outer Dimensions HxWxD, mm | 870x555x555 |
| Inner Dimensions HxWxD, mm | 465x390x390 |
| Weight Gross/Net, kg | 56 / 52 |
| Material inner cabinet | Painted Stee |
| Material outer cabinet | Painted Stee |
| Packaging weight, kg | |
| Packaging dimensions HxWxD, mm | 1055x675x625 |
| Insulation thickness, mm | 80 |
| Insulation type | Polyurethane with Cyclopentane |
| Mobility / Castors | No |
| Refrigerant, Type | Nature F |
| Number of compressors | |
| Internal Air Distribution | Statio |
| Number of Probes | |
| CONTROLLER | |
| Controller | XR30C> |
| Controller language | No language - only 3 digits |
| USB Connection | No. |
| Logging | No |
| Temperature Graph | No |
| High/Low Temp. Alarm | Yes |
| Open Door Alarm | No |
| Probe Failure Alarm | Yes |
| Power Failure Alarm | No |
| STORAGE | |
| Volume, Gross/net, L | 7/2 |
| Baskets | |
| Basket material | Steel coated with plastic powder |
| Inner lids | Nc |
| FEATURES | |
| Lock | Yes |
| LED Light | No |
| Battery Backup for Controller, 24h | No |
| Porthole | Yes - Ø 12 mm |
| Dry Contact | No |
| Door | Solic |
| Door Reversibility | N/A |

SOLUTIONS



BIOMEDICAL

The VT low temperature freezers creates the possibility to maintain temperatures as low as -60°C. Supreme stability, reliability, user-friendliness and ease of cleaning make these freezers an ideal solution for laboratories and hospitals.

| Frequency | Hz | 50Hz |
|--|----------|------------|
| Max Ambient | °C | 30°C |
| Max Humidity | % rh | 55% |
| PERFORMANCE | | |
| All data in RT20°C | | |
| Temperature Range | °C | -40 to -60 |
| Uniformity in performance - difference +/- from Avg set point | °C | - |
| Pull down time (from +25° to -40°) | Minutes | 60 |
| Hold over time (from factory SP to -25, -40 and -60) Empty, In 25° | Minutes | 80 |
| Refrigerant | | Nature R |
| Number of probes | pcs | 1 |
| Defrost | y/n | No |
| Internal air distribution | | Static |
| Number of compressors | pcs | 1 |
| Safety thermostat | y/n | No |
| Energy 24 hours, in 25° | kWh/24h | 2,24 |
| Energy year | kWh/year | 818 |

2



VT 75 -60°C LOW TEMPERATURE FREEZER

| Height, mm | 895 |
|------------------|-----|
| Width, mm | 550 |
| Depth, mm | 550 |
| Incl. handle, mm | 645 |

All measurements in mm.



